Electro Polished Stainless Steel Chamber (D shaped Box)

4" Diameter View Port on Front Door with Manual Shutter

Single/Multiple Magnetron Gun Assemblies with Various Target Sizes

Manual Shutter Assembly for All Sources

DC / RF Sputtering Deposition

Turbo Molecular Vacuum Pumping System with Matching Dual Stage Rotary Vane Pump


Quartz Crystal Thickness Sensor with Deposition Controller

Mass Flow Controller with Digital Readout

Full Range Vacuum Gauge with Digital Display & Readout

Easy Target Mounting

PLC Controlled System



Water Chiller

Substrate Pre-Cleaning

Load Lock Subsystem with Sample Transfer Adaptability

Multi-Axis Substrate Stage Rotation

Quartz Lamp Heater (from 300° C up to 800° C)

Planetary Substrate Stage

PC Controlled System

Water Cooled Substrate Stage

Adjustable Size & Height Available for Substrate Stage SQC 310 Film Thickness Monitor &

Deposition Controller

Motorized Shutter Assembly

Dry Scroll Pump

Cryo Pumping System

Cold Cathode or Hot Cathode Ionization Gauge

Additional Spare Flanges for Future Upgrades



TORR INTERNATIONAL, INC. THIN FILM & NANOTECHNOLOGY

Email: torr@torr.com www.torr.com (845) 565-4027 (888) MAC-TORR Fax: (845) 561-7731

1 Torr = 1 mm Ha